

INO

QuickPOZ

Optomechanical Mounts and Breadboards 2022-2023 Catalog

EC19500 / INO-PLT0201-0037 version 5.0



Table of Content

QuickPOZ, INO's family of robust optomechanical mounts for accurate positioning	3
How does it work?	3
Performance specifications	5
What is included in INO's QuickPOZ family?	6
What can be done with these mounts?	9
Need a custom configuration?	9
General considerations.....	10
Tooling & miscellaneous	11
Reference ball, 6.35mm	12
Threaded ring	13
Split ring	14
XY adjustment tooling, D05 & D1	16
Breadboard locating union	17
Breadboard clamp tool	18
Mirror mount leaf spring assembly	19
Male-female QC thread adaptor	21
Threaded iris	22
Optical spacer ring	23
Optical spacer, special order	24
Mounting Plates	25
Breadboard	26
Breadboard disc, gimbal	27
Gimbal adaptor for mirror mount and fixed mount series	28
Mirror Mounts.....	29
Fixed mirror mount, horizontal	30
Fixed mirror mount, vertical-bottom	32
Adjustable mirror mount, horizontal	34
Generic Mounts.....	36
XY adjustable mount	37
Through-hole XY adjustable mount	38
XYZ adjustable mount	39

Fixed mount with flange	40
Fixed mount	41
Tubes.....	42
Lens tube, male-female	43
Lens tube, male-female, adapted for MA_TH_	44
Lens tube with axial adjustment	45
Rotation Mounts.....	46
Rotation mount	47
Translation Mounts.....	48
Translation mount, 1.5mm, 62.5MM	49

QuickPOZ, INO's family of robust optomechanical mounts for accurate positioning

The QuickPOZ optical mounts and breadboard line-up is addressing the need for robust optomechanical prototyping. Using these self-positioning optomechanical mounts is a cost-effective way of rapidly assembling prototypes that will remain aligned even under severe operating or transportation conditions. All mounts have been designed to be operated under a typical transport vehicle vibration environment up to 500 Hz while keeping their pointing stability under $\pm 50 \mu\text{rad}$.

The nominal positioning of all optics on the breadboard is easy with INO's QuickPOZ, since each mount can be located with high position repeatability using removeable reference balls on INO's special breadboard. These mounts integrate the patented QuickCTR-edge technology (US Patent 11,327,332) allowing the positioning of all optical components within $\pm 50 \mu\text{m}$ @ 2 sigma from any mount assembled on the same breadboard with respect to nominal optical axis.

There are three standard heights for the mount optical axis with respect to the breadboard: 25.4 mm, 31.8 mm, and 38.1 mm.

How does it work?

INO's QuickPOZ is the combination of robust optomechanical mounts and a patented technology to accurately and rapidly position them on a breadboard.

These optomechanical mounts are the fruit of more than 25 years of development of prototypes used in demanding environments. The mounts included in the QuickPOZ line-up are found in surveillance planes, severe industrial environments, and bioscience labs – to name a few. High accuracy threaded reference balls are temporarily installed on INO's proprietary breadboard to locate each mount with respect to one another. Breadboard threaded holes are accurately manufactured allowing precise positioning of all mounts. To complement the component positioning, each mount has been designed with tight tolerance reference features.




Figure 1 QuickPOZ breadboard reference ball positioning method

The lenses are autocentered into the mounts using INO's patented QuickCTR-edge technology which is based on the geometrical relationship between the lens chamfer and the threaded ring contact seat radius.

If needed, mounts from other vendors can be fitted on QuickPOZ breadboards.




Figure 2 QuickCTR-edge technology principle

The technology is applicable to multiple mounting configurations like: convex, concave, and plano optical surfaces; optical subassemblies; and tube stacks.




Figure 3 A few configurations using QuickCTR-edge technology

Performance specifications





Feature	Performance
Optomechanical mount compatibility	Compatible with standard threads 0.535"-40, 1.035"-40, 2.035"-40 and RMS commercial accessories (Thorlabs, Edmund, Newport, ...). When a commercial accessory is used, centering and positioning performances are not met. In some cases, their robustness may also not be met.
Operating temperature	-40°C to +50°C
Storage temperature	-46°C to +63°C MIL-STD-810H Method 501.7 Procedure I, minimum 7 cycles (25°C to 63°C) with 2-hour plateau, 3°C/min MIL-STD-810H Method 502.7 Procedure I, 1 cycle (25°C to -46°C), 24-hour plateau, 3°C/min
Optical axis positioning accuracy	Statistical RSS accuracy of ± 0.05 mm between optical axis of any optomechanical mount with respect to the nominal optical axis of the breadboard assembly. For breadboard dimensions up to 200 x 400 mm or 300 x 300 mm.
Mirror surface deformations induced by mount	MMA, MMH, and MMV mirror mount series induced deformations are $\leq \lambda/10$ PV @ 633 nm over their clear aperture. Measured on $\varnothing 25.4$ mm x 6.13 mm and $\varnothing 50.8$ x 12 mm mirror substrates mounted with a stack of 3x3 blades with a maximal deflection of 0.5 mm. ***Do not exceed 0.5 mm of blade deflection to avoid creeping***
Dimensional pointing stability over operating temperature range	$\leq \pm 50$ μrad (mechanical angle) Optics below $\varnothing 25$ mm may exceed the ± 50 μ rad pointing stability due to their small size and small mounting seat.
Mounting repeatability	$\leq \pm 0.015$ mm in positioning
Shipping vibrations impact on angular positioning	Without shipping packaging: $\leq \pm 50$ μrad (mechanical angle). MIL-STD-810H Method 514.8 C.II category 4, unknown orientation, random vibration, 20 -500Hz, 1.17 Grms. Within a typical cardboard shipping packaging: $\leq \pm 50$ μrad (mechanical angle). MIL-STD-810H Method 514.8 E-1 category 24, all axis orientation, random vibration, 20 -2000Hz, 7.7 Grms, 1 hour/axis.
Shock survival	30 G minimum, without shipping packaging. 30 G corresponds to the limit of the most sensitive components, which are $\varnothing 50.8$ mm mirror mounts.
Stress relief	If necessary, it is possible to improve dimensional stability performances by conditioning the assembly to thermal stress relief cycles. Thermal stress relief cycles are application specific and may be available upon request.


What is included in INO's QuickPOZ family?

The whole family covers nearly 150 mounts of several sizes (QC05, QC1, QC30, QC2), available in three optical axis heights (25.4 mm, 31.8 mm, and 38.1 mm), and in four different mirror diameters. Adjustable mounts are also available for transversal (normal to optical axis), axial, tip/tilt, and clocking positioning.

QuickPOZ mount and accessory overview

	Product name	Description
	Threaded rings	Autocentered threaded rings for optical components, available in QC05, QC1, QC30, and QC2 thread sizes.
	Adjustment tools	Several removable tool designed to be used with the QuickPOZ mounts for nominal or precise alignment..
	Threaded iris	Autocentered Ø1 mm iris, available in QC05 & QC1 thread sizes. Used for alignment purposes.
	Fixed mirror mounts, horizontal	Low distortion fixed mirror mounts for horizontal beam folding, available in Ø12.7 mm, Ø25.4 mm, Ø38.1 mm, and Ø50.8 mm sizes.
	Fixed mirror mounts, vertical-bottom	Low distortion fixed mirror mounts for downwards beam folding, available in Ø12.7 mm, Ø25.4 mm, Ø38.1 mm, and Ø50.8 mm sizes.

	Adjustable mirror mounts	<p>Low distortion adjustable mirror mounts with $\pm 2^\circ$ tip-tilt for horizontal beam folding, available in $\varnothing 12.7$ mm, $\varnothing 25.4$ mm, $\varnothing 38.1$ mm, and $\varnothing 50.8$ mm sizes.</p> <p>Left-hand and right-hand versions are available.</p>
	XY adjustable mounts	X-Y adjustable mounts (± 1 mm) with QC05 & QC1 threads and axial support, for submicron positioning with removable XY adjustment tool.
	XYZ adjustable mounts	X-Y-Z adjustable mounts with QC1 threads, for submicron positioning with removable XY adjustment tool (± 1 mm), and 7mm travel in Z.
	Fixed mounts	Fixed mounts with and without flange, available in autocentered thread sizes QC05, QC1, QC30, and QC2, and microscope size 0.8"-36 (RMS).
	Lens tubes	<p>Standard autocentered tubes of several lengths, stackable, available sizes.</p> <p>Some can be used with through-hole XY adjustable mounts.</p> <p>Tube thread adaptors and autocentered axially adjustable tubes are also available.</p>
	Rotation mounts	Compact rotating mounts with autocentered thread sizes QC05 & QC1.

	Translation mounts	Ultra-stable and high accuracy translation stage (± 1.5 mm) to be combined with adjustable and fixed mounts.
---	--------------------	---

Correspondence between INO's autocentered threads and industry standards

QuickPOZ thread	Industry standard equivalent
QC05	0.535"-40
QC1	1.035"-40
QC30	<i>This thread is not compatible with industry standard.</i>
QC2	2.035"-40
RMS	RMS

Many mounting breadboards, plates, and associated accessories are also available.




Figure 4 QuickPOZ mounting plates overview

What can be done with these mounts?

Any industrial optical applications where robustness and precise positioning are a concern may benefit from QuickPOZ. There are infinite ways of using these mounts, whether for an industrial laser source, a bulk fibre optical device, a spectrometer, a medical optical device, an illumination system, or an objective lens.



Figure 5 First example of configuration




Figure 6 Second example of configuration

Need a custom configuration?

INO is offering consultation services to build up your own configuration for optimum performances. INO can also fully customize your need, from a custom breadboard up to a full turnkey solution. Our specialists in optical and optomechanical design can assist you during the whole design development process and even for your production series.

General considerations

All QuickPOZ mounts are designed to be attached using M4x0.7 screws which are 14 mm or 20 mm in length, depending of the mount type. In all cases, a washer Ø9 mm x 0.8 mm thick (McMaster #93475A230) must be used with the deburred side downwards to avoid damaging the mount.

QuickPOZ specifications are guaranteed only if screws have the proper tightening torque:

Screw Size Diameter (mm)	Thread Pitch (mm)	Torque			
		(N-m)	(ozf-in)	(lbf-in)	(lbf-ft)
1.6	0.35	0.12	17		
2	0.4	0.25	35		
2.5	0.45	0.51	72	4	
3	0.5	0.91	128	8	
4	0.7	2.11		19	
5	0.8	4.26		38	
6	1	7.24		64	
8	1	18.82			14
10	1.25	36.72			27
14	1.5	104.58			77

All QuickPOZ threaded rings are compatible with Thorlabs spanner wrenches for SM05, SM1, SM2, and SM30 formats.


Tooling & miscellaneous

Reference ball, 6.35mm

Description	Positioning ball assembly
Adjustments	n/a
Required tool	2mm Allen wrench, or fingers
Product notes	Only lightly tighten by hand; only use hex key for accessibility.

Part no.	Unit Price (\$)
TLBAL1	TBA

Drawing TLBAL1



Threaded ring

Description	Threaded ring for autocentered optical components
Adjustments	n/a
Required tool	Compatible with Thorlabs spanner wrench for ring series SM05RR, SM1RR, SM2RR, SM30RR
Product notes	Compatible with QuickPOZ QC_ and commercial tube series.

Part no.	Unit Price (\$)
RQC05	TBA
RQC1	TBA
RQC30	TBA
RQC2	TBA


Drawing RQC_

<p>4 3 2 1</p> <p>PROPRIETARY NOTICE/CONFIDENTIALITY: THIS DOCUMENT AND THE INFORMATION DISCLOSED HEREIN ARE CONFIDENTIAL AND THE PROPERTY OF INSTITUT NATIONAL D'OPTIQUE (NATIONAL OPTICS INSTITUTE, HEREINAFTER "INO"). THEY SHALL NOT BE REPRODUCED OR DISCLOSED IN WHOLE OR IN PART WITHOUT THE EXPRESS WRITTEN AUTHORIZATION OF INO. INO RESERVES ALL PATENT, COPYRIGHT AND OTHER PROPRIETARY RIGHTS TO THIS DOCUMENT, INCLUDING ALL DESIGN, MANUFACTURING, REPRODUCTION, USE AND SALES RIGHTS THEREON, REPRODUCTION, USE AND SALES RIGHTS THEREON, EXCEPT TO THE EXTENT SAID RIGHTS ARE EXPRESSLY GRANTED TO OTHERS.</p> <p>NOTES: 1 - FOR REFERENCE ONLY, NOT FOR MANUFACTURING PURPOSES</p>				<p>REVISIONS</p> <table border="1"> <tr> <th>REV.</th> <th>DESCRIPTION</th> <th>DATE (YYYY-MM-DD)</th> <th>DRAWN</th> </tr> <tr> <td> </td> <td> </td> <td> </td> <td> </td> </tr> </table>				REV.	DESCRIPTION	DATE (YYYY-MM-DD)	DRAWN																															
REV.	DESCRIPTION	DATE (YYYY-MM-DD)	DRAWN																																							
<p>THK</p> <p>WIDTH</p> <p>Ø ID</p> <p>THREAD</p>																																										
<table border="1"> <thead> <tr> <th>PRODUCT #</th> <th>NAME</th> <th>WIDTH</th> <th>THREAD</th> <th>ID</th> <th>THK</th> <th>MASS (KG)</th> </tr> </thead> <tbody> <tr> <td>RQC05</td> <td>THREADED RING QC05</td> <td>1.191</td> <td>0.535"-40 CLASS 3A</td> <td>11.00</td> <td>2.480</td> <td>0.0003</td> </tr> <tr> <td>RQC1</td> <td>THREADED RING QC1</td> <td>0.794</td> <td>1.035"-40 CLASS 3A</td> <td>22.90</td> <td>2.781</td> <td>0.0008</td> </tr> <tr> <td>RQC30</td> <td>THREADED RING QC30</td> <td>0.794</td> <td>M31 X 0.5 - 4h</td> <td>27.94</td> <td>2.781</td> <td>0.0009</td> </tr> <tr> <td>RQC2</td> <td>THREADED RING QC2</td> <td>0.794</td> <td>2.035"-40 CLASS 3A</td> <td>48.26</td> <td>3.281</td> <td>0.0020</td> </tr> </tbody> </table>				PRODUCT #	NAME	WIDTH	THREAD	ID	THK	MASS (KG)	RQC05	THREADED RING QC05	1.191	0.535"-40 CLASS 3A	11.00	2.480	0.0003	RQC1	THREADED RING QC1	0.794	1.035"-40 CLASS 3A	22.90	2.781	0.0008	RQC30	THREADED RING QC30	0.794	M31 X 0.5 - 4h	27.94	2.781	0.0009	RQC2	THREADED RING QC2	0.794	2.035"-40 CLASS 3A	48.26	3.281	0.0020	<p>UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN MILLIMETERS REMOVE BURRS AND SHARP EDGES TOLERANCES AFTER SURFACE TREATMENT 0.8/ SHOWN: ±0.1, ∠±14° DRAWING IN ACCORDANCE WITH ASME Y14.5 2009</p> <p>CONTRACT NUMBER/INITIAL PROJECT 189810</p> <p>APPROVALS DESIGNED F LAMONTAGNE DATE (YYYY-MM-DD) 2018-11-26</p> <p>DRAWN J. RÉGNIER DATE 2020-03-02</p> <p>EXAM/CHECK <i>PR</i> DATE 2020-03-30</p> <p>VÉRIF./VERIF.: <i>MD</i> DATE 2020-03-30</p> <p>MATERIAL ALUMINUM</p> <p>SURFACE TREATMENT BLACK ANODIZING</p> <p>INO ALL RIGHTS RESERVED © NATIONAL OPTICS INSTITUTE 2740, EINSTEIN STREET QUEBEC CITY, QUEBEC CANADA, G1P 4S4 T (418) 657-7006 F (418) 657-7009 www.ino.ca</p> <p>TITLE THREADED RING</p> <p>FORMAT A DAI (CAGE CODE) 3AT79</p> <p>DRAWING NUMBER RQC_</p> <p>REVISION A</p> <p>SCALE 2:1</p> <p>MASS (kg)</p> <p>PAGE 1/1</p>			
PRODUCT #	NAME	WIDTH	THREAD	ID	THK	MASS (KG)																																				
RQC05	THREADED RING QC05	1.191	0.535"-40 CLASS 3A	11.00	2.480	0.0003																																				
RQC1	THREADED RING QC1	0.794	1.035"-40 CLASS 3A	22.90	2.781	0.0008																																				
RQC30	THREADED RING QC30	0.794	M31 X 0.5 - 4h	27.94	2.781	0.0009																																				
RQC2	THREADED RING QC2	0.794	2.035"-40 CLASS 3A	48.26	3.281	0.0020																																				

Split ring

Description	Complementary ring for lens with small convex radius of curvature
Adjustments	n/a
Required tool	n/a
Product notes	<p>This split ring needs to be added between the lens and the QuickPOZ Threaded Ring when the radius of curvature of a convex lens is smaller than the threaded ring mounting radius.</p> <p>Compatible with QuickPOZ QC_ and commercial tube series.</p>

Part no.	Unit Price (\$)
SRD05	TBA
SRD1	TBA
SRD30	TBA
SRD2	TBA



This split ring is required for convex lens radius of curvature smaller than:

Threaded ring P/N	Max. radius of curvature of convex lens (mm)
RQC05	11.5
RQC1	23.0
RQC30	27.5
RQC2	47.0

Drawing SRD_ (...see next page)


4	3	2	1																														
<small>PROPRIETARY NOTICE/CONFIDENTIALITY: THIS DOCUMENT AND THE INFORMATION DISCLOSED HEREIN ARE CONFIDENTIAL AND THE PROPERTY OF INSTITUT NATIONAL D'OPTIQUE NATIONAL OPTICS INSTITUTE (HEREAFTER "INO"). THEY SHALL NOT BE REPRODUCED OR DISCLOSED IN WHOLE OR IN PART WITHOUT THE EXPRESS WRITTEN AUTHORIZATION OF INO. INO RESERVES ALL PATENT, COPYRIGHT AND OTHER PROPRIETARY RIGHTS TO THIS DOCUMENT, INCLUDING ALL DESIGN, MANUFACTURING, REPRODUCTION, USE AND SALES RIGHTS THEREON, REPRODUCTION, USE AND SALES RIGHTS THEREON, EXCEPT TO THE EXTENT SAID RIGHTS ARE EXPRESSLY GRANTED TO OTHERS.</small>		REVISIONS																															
NOTES: 1 - FOR REFERENCE ONLY, NOT FOR MANUFACTURING PURPOSES		REV.	DESCRIPTION																														
		DATE (YYYY-MM-DD)	DRAWN																														
<div style="display: flex; justify-content: space-around;"> <div style="text-align: center;"> <p>Ø OD Ø ID</p> </div> <div style="text-align: center;"> <p>THK</p> </div> <div style="text-align: center;"> <p>SRD05</p> <p>SEE DETAIL A</p> </div> </div>	<div style="display: flex; justify-content: space-around;"> <div style="text-align: center;"> <p>FOR SRD05 ONLY : INSTALLATION NOT REVERSIBLE. THIS SIDE MUST BE IN CONTACT WITH RETAINING RING RQC05.</p> <p>DETAIL A SCALE 5:1</p> </div> <div style="text-align: center;"> </div> </div>																																
<table border="1" style="width: 100%; border-collapse: collapse; font-size: 0.8em;"> <thead> <tr> <th>PRODUCT #</th> <th>NAME</th> <th>ID</th> <th>OD</th> <th>THK</th> <th>MASS (KG)</th> </tr> </thead> <tbody> <tr> <td>SRD05</td> <td>SPLIT RING, D05</td> <td>11.0</td> <td>12.98</td> <td>2.50</td> <td>0.0002</td> </tr> <tr> <td>SRD1</td> <td>SPLIT RING, D1</td> <td>22.9</td> <td>25.68</td> <td>3.00</td> <td>0.0008</td> </tr> <tr> <td>SRD30</td> <td>SPLIT RING, D30</td> <td>27.9</td> <td>30.56</td> <td>3.25</td> <td>0.0010</td> </tr> <tr> <td>SRD2</td> <td>SPLIT RING, D2</td> <td>48.3</td> <td>51.08</td> <td>4.00</td> <td>0.0022</td> </tr> </tbody> </table>				PRODUCT #	NAME	ID	OD	THK	MASS (KG)	SRD05	SPLIT RING, D05	11.0	12.98	2.50	0.0002	SRD1	SPLIT RING, D1	22.9	25.68	3.00	0.0008	SRD30	SPLIT RING, D30	27.9	30.56	3.25	0.0010	SRD2	SPLIT RING, D2	48.3	51.08	4.00	0.0022
PRODUCT #	NAME	ID	OD	THK	MASS (KG)																												
SRD05	SPLIT RING, D05	11.0	12.98	2.50	0.0002																												
SRD1	SPLIT RING, D1	22.9	25.68	3.00	0.0008																												
SRD30	SPLIT RING, D30	27.9	30.56	3.25	0.0010																												
SRD2	SPLIT RING, D2	48.3	51.08	4.00	0.0022																												
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN MILLIMETERS REMOVE BURRS AND SHARP EDGES TOLERANCES AFTER SURFACE TREATMENT 0.8 SHOWN: ±0.1, < ±1/4" DRAWING IN ACCORDANCE WITH ASME Y14.5 2009		CONTRACT NUMBER/INITIAL PROJECT 189810																															
MATERIAL ALUMINUM		INO NATIONAL OPTICS INSTITUTE 2740, EINSTEIN STREET QUEBEC CITY, QUEBEC CANADA, G1P 4S4 T (418) 657-7006 F (418) 657-7009 www.ino.ca																															
SURFACE TREATMENT BLACK ANODIZING		TITLE SPLIT RING																															
APPROVALS DESIGNED F. LAMONTAGNE DRAWN J. RÉGNIER VERIFIED APPROVED		DATE (YYYY-MM-DD) 2019-10-10 2020-03-30																															
FORMAT DAI (CAGE CODE) A 3AT79		DRAWING NUMBER SRD																															
SCALE 2:1		MASS (kg) PAGE 1/1																															

XY adjustment tooling, D05 & D1

Description	Transverse X-Y manipulator, removable with differential screws; fits with QuickPOZ mounts MA_QC05, MA_QC1, MA_TH05, MA_TH1, and MA_TAQC1.
Adjustments	$\pm 1.5\text{mm}$; coarse $318\mu\text{m/rev.}$ and fine $25\mu\text{m/rev.}$
Required tool	2mm Allen wrench
Product notes	<p>Can be used either with the adjuster located at right or at left.</p> <p>***Warning*** Do not forget to detent the spring plunger before adjustment (small set screws located in front of the mount).</p> <p>***Warning*** Do not forget to retract the spring plunger before removing the tool from the mount.</p>

Part no.	Unit Price (\$)
TLXY1	TBA

Drawing TLXY1




Breadboard locating union

Description	Tool to join 2 mounting plates together
Adjustments	n/a
Required tool	n/a
Product notes	User instructions available upon request

Part no.	Unit Price (\$)
TLBBU	TBA

Drawing TLBBU



Breadboard clamp tool

Description	Clamp tool to fix breadboard gimbals, or to fix mounts and other breadboards (6.4mm, 7.6mm, or 12.7mm thick).	Part no.	Unit Price (\$)
Adjustments	n/a	TLCP064	TBA
Required tool	n/a	TLCP076	TBA
Product notes		TLCP127	TBA

Drawing TLCP_

<p>4 3 2 1</p> <p>PROPRIETARY NOTICE/CONFIDENTIALITY: THIS DOCUMENT AND THE INFORMATION DISCLOSED HEREIN ARE CONFIDENTIAL AND THE PROPERTY OF INSTITUT NATIONAL D'OPTIQUE NATIONAL OPTICS INSTITUTE (HEREAFTER INO). THEY SHALL NOT BE REPRODUCED OR DISCLOSED IN WHOLE OR IN PART WITHOUT THE EXPRESS WRITTEN AUTHORIZATION OF INO. INO RESERVES ALL PATENT, COPYRIGHT AND OTHER PROPRIETARY RIGHTS TO THIS DOCUMENT, INCLUDING ALL DESIGN, MANUFACTURING, REPRODUCTION, USE AND SALES RIGHTS THEREIN. REPRODUCTION, USE AND SALES RIGHTS THEREIN, EXCEPT TO THE EXTENT SAID RIGHTS ARE EXPRESSLY GRANTED TO OTHERS.</p> <p>NOTES: 1 - FOR REFERENCE ONLY, NOT FOR MANUFACTURING PURPOSES</p>				<p>REV. DESCRIPTION DATE (YYYY-MM-DD) DRAWN</p>																							
<p>Ø 4.8</p> <p>17.98</p> <p>9</p>																											
<p>HEIGHT</p> <p>CLAMP SIZE</p>				<table border="1"> <thead> <tr> <th>PRODUCT #</th> <th>NAME</th> <th>HEIGHT</th> <th>CLAMP SIZE</th> <th>MASS (KG)</th> </tr> </thead> <tbody> <tr> <td>TLCP064</td> <td>CLAMP, 6.35MM</td> <td>9.600</td> <td>6.350</td> <td>0.002</td> </tr> <tr> <td>TLCP076</td> <td>CLAMP, 7.6MM</td> <td>10.850</td> <td>7.600</td> <td>0.002</td> </tr> <tr> <td>TLCP127</td> <td>CLAMP, 12.7MM</td> <td>15.950</td> <td>12.700</td> <td>0.003</td> </tr> </tbody> </table>				PRODUCT #	NAME	HEIGHT	CLAMP SIZE	MASS (KG)	TLCP064	CLAMP, 6.35MM	9.600	6.350	0.002	TLCP076	CLAMP, 7.6MM	10.850	7.600	0.002	TLCP127	CLAMP, 12.7MM	15.950	12.700	0.003
PRODUCT #	NAME	HEIGHT	CLAMP SIZE	MASS (KG)																							
TLCP064	CLAMP, 6.35MM	9.600	6.350	0.002																							
TLCP076	CLAMP, 7.6MM	10.850	7.600	0.002																							
TLCP127	CLAMP, 12.7MM	15.950	12.700	0.003																							
<p>UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN MILLIMETERS REMOVE BURRS AND SHARP EDGES TOLERANCES AFTER SURFACE TREATMENT 0.8 SHOWN: ±0.1, ∠±14° DRAWING IN ACCORDANCE WITH ASME Y14.5 2009</p>				<p>CONTRACT NUMBER/INITIAL PROJECT 189810</p>																							
<p>MATERIAL ALUMINUM</p>				<p>APPROVALS DATE (YYYY-MM-DD)</p> <p>DESIGNED S. PARADIS 2019-06-10</p>																							
<p>SURFACE TREATMENT</p>				<p>DRAWN S. PARADIS 2019-10-07</p>																							
<p>BLACK ANODIZING</p>				<p>EXAM/CHECK <i>PR</i> 2020-03-17</p>																							
<p>VÉRIF./VERIF.: <i>MD</i> 2020-03-17</p>				<p>FORMAT DAI (CAGE CODE) DRAWING NUMBER</p> <p>A 3AT79 TLCP_</p>																							
<p>SCALE 2:1</p>				<p>MASS (kg) PAGE 1/1</p>																							
<p>4 3 2 1</p>				<p>10 TLCP_</p>																							

Mirror mount leaf spring assembly

Description	3 leaf springs with 1 captive screw
Adjustments	n/a
Required tool	2mm Allen wrench
Product notes	<p>This stack of 3 leaf springs gives the following force per mounting point:</p> <ul style="list-style-type: none"> • 5.9 ±0.5 N @ 0.5mm deflection • 3.9 ±0.5 N @ 0.4mm deflection • 2.9 ±0.2 N @ 0.3mm deflection • 2.0 ±0.2 N @ 0.2mm deflection <p>***Warning*** These leaf springs have been designed to be used in stack of 3 with a maximal deflection of 0.5mm at the tip.</p>

Part no.	Unit Price (\$)
MMLSA	TBA

Drawing MMLSA


4	3	2	1																								
<p><small>PROPRIETARY NOTICE/CONFIDENTIALITY: THIS DOCUMENT AND THE INFORMATION DISCLOSED HEREIN ARE CONFIDENTIAL AND THE PROPERTY OF INSTITUT NATIONAL D'OPTIQUE (NATIONAL OPTICS INSTITUTE, HEREINAFTER "INO"). THEY SHALL NOT BE REPRODUCED OR DISCLOSED IN WHOLE OR IN PART WITHOUT THE EXPRESS WRITTEN AUTHORIZATION OF INO. INO RESERVES ALL PATENT, COPYRIGHT AND OTHER PROPRIETARY RIGHTS TO THIS DOCUMENT, INCLUDING ALL DESIGN, MANUFACTURING, REPRODUCTION, USE AND SALES RIGHTS THEREON. REPRODUCTION, USE AND SALES RIGHTS THEREON, EXCEPT TO THE EXTENT SAID RIGHTS ARE EXPRESSLY GRANTED TO OTHERS.</small></p> <p>NOTES: 1 - FOR REFERENCE ONLY, NOT FOR MANUFACTURING PURPOSES</p>		<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th colspan="4" style="text-align: center;">REVISIONS</th> </tr> <tr> <th style="width: 10%;">REV.</th> <th style="width: 60%;">DESCRIPTION</th> <th style="width: 20%;">DATE (YYYY-MM-DD)</th> <th style="width: 10%;">DRAWN</th> </tr> <tr> <td> </td> <td> </td> <td> </td> <td> </td> </tr> </table>		REVISIONS				REV.	DESCRIPTION	DATE (YYYY-MM-DD)	DRAWN																
REVISIONS																											
REV.	DESCRIPTION	DATE (YYYY-MM-DD)	DRAWN																								
<p>UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN MILLIMETERS REMOVE BURRS AND SHARP EDGES TOLERANCES AFTER SURFACE TREATMENT 0.8 SHOWN: ± 0.1, $\angle \pm 1/4^\circ$ DRAWING IN ACCORDANCE WITH ASME Y14.5 2009</p> <p>MATERIAL: STEEL</p> <p>SURFACE TREATMENT: N/A</p>		<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="2" style="text-align: center;">CONTRACT NUMBER/INITIAL PROJECT 189810</td> <td colspan="2" style="text-align: center;"> <p><small>NATIONAL OPTICS INSTITUTE 2740, EINSTEIN STREET QUEBEC CITY, QUEBEC CANADA, G1P 4S4</small></p> <p><small>T (418) 657-7006 F (418) 657-7009 www.ino.ca</small></p> </td> </tr> <tr> <td colspan="2" style="text-align: center;">APPROVALS</td> <td colspan="2" style="text-align: center;">DATE (YYYY-MM-DD)</td> </tr> <tr> <td>DESIGNED</td> <td>B. LEDUC</td> <td colspan="2">2019-10-01</td> </tr> <tr> <td>DRAWN</td> <td>S. PARADIS</td> <td colspan="2">2020-03-24</td> </tr> <tr> <td>EXAM./CHECK</td> <td><i>SP</i></td> <td colspan="2">2020-03-30</td> </tr> <tr> <td>VÉRIF./VÉRIF.</td> <td><i>MD</i></td> <td colspan="2">2020-03-30</td> </tr> </table>		CONTRACT NUMBER/INITIAL PROJECT 189810		<p><small>NATIONAL OPTICS INSTITUTE 2740, EINSTEIN STREET QUEBEC CITY, QUEBEC CANADA, G1P 4S4</small></p> <p><small>T (418) 657-7006 F (418) 657-7009 www.ino.ca</small></p>		APPROVALS		DATE (YYYY-MM-DD)		DESIGNED	B. LEDUC	2019-10-01		DRAWN	S. PARADIS	2020-03-24		EXAM./CHECK	<i>SP</i>	2020-03-30		VÉRIF./VÉRIF.	<i>MD</i>	2020-03-30	
CONTRACT NUMBER/INITIAL PROJECT 189810		<p><small>NATIONAL OPTICS INSTITUTE 2740, EINSTEIN STREET QUEBEC CITY, QUEBEC CANADA, G1P 4S4</small></p> <p><small>T (418) 657-7006 F (418) 657-7009 www.ino.ca</small></p>																									
APPROVALS		DATE (YYYY-MM-DD)																									
DESIGNED	B. LEDUC	2019-10-01																									
DRAWN	S. PARADIS	2020-03-24																									
EXAM./CHECK	<i>SP</i>	2020-03-30																									
VÉRIF./VÉRIF.	<i>MD</i>	2020-03-30																									
<p>4</p>		<p>1</p>																									

Male-female QC thread adaptor

Description	Male-female thread adaptor
Adjustments	n/a
Required tool	n/a
Product notes	

Part no.	Unit Price (\$)
AMFQC05QC1	TBA
AMFQC05QC30	TBA
AMFQC1QC05	TBA
AMFQC1QC30	TBA
AMFQC1QC2	TBA
AMFQC30QC05	TBA
AMFQC30QC1	TBA
AMFQC2QC1	TBA

Drawing AMFQC_QC_




Threaded iris

Description	Autocentered threaded iris, QC1 or QC05, with a 1mm hole
Adjustments	n/a
Required tool	n/a
Product notes	These iris are used for alignment purposes

Part no.	Unit Price (\$)
TI100QC05	TBA
TI100QC1	TBA

Drawing TI_QC_



Optical spacer ring

Description Shims to be used with QuickPOZ mirror mount series, to fill in the thickness gap obtained with some mirrors, filters, or dichroics

Adjustments n/a

Required tool n/a

Product notes

Part no.	Unit Price (\$)	Part no.	Unit Price (\$)
OST05D05	TBA	OST05D15	TBA
OST1D05	TBA	OST1D15	TBA
OST2D05	TBA	OST2D15	TBA
OST3D05	TBA	OST3D15	TBA
OST5D05	TBA	OST5D15	TBA
OST05D1	TBA	OST05D2	TBA
OST1D1	TBA	OST1D2	TBA
OST2D1	TBA	OST2D2	TBA
OST3D1	TBA	OST3D2	TBA
OST5D1	TBA	OST5D2	TBA

Drawing OST_D_


4		3		2		1																																																																																																																															
<p>PROPRIETARY NOTICE/CONFIDENTIALITY: THIS DOCUMENT AND THE INFORMATION DISCLOSED HEREIN ARE CONFIDENTIAL AND THE PROPERTY OF INSTITUT NATIONAL D'OPTIQUE (INO). NO PART OF THIS DOCUMENT SHALL BE REPRODUCED OR DISCLOSED IN WHOLE OR IN PART WITHOUT THE EXPRESS WRITTEN AUTHORIZATION OF INO. INO RESERVES ALL PATENT, COPYRIGHT AND OTHER PROPRIETARY RIGHTS TO THIS DOCUMENT, INCLUDING ALL DESIGN, MANUFACTURING, REPRODUCTION, USE AND SALES RIGHTS THEREIN. REPRODUCTION, USE AND SALES RIGHTS THEREIN, EXCEPT TO THE EXTENT SAID RIGHTS ARE EXPRESSLY GRANTED TO OTHERS.</p> <p>NOTES: 1 - FOR REFERENCE ONLY, NOT FOR MANUFACTURING PURPOSES</p>				<p>REVISIONS</p> <table border="1"> <tr> <th>REV.</th> <th>DESCRIPTION</th> <th>DATE (YYYY-MM-DD)</th> <th>DRAWN</th> </tr> <tr> <td> </td> <td> </td> <td> </td> <td> </td> </tr> </table>				REV.	DESCRIPTION	DATE (YYYY-MM-DD)	DRAWN																																																																																																																										
REV.	DESCRIPTION	DATE (YYYY-MM-DD)	DRAWN																																																																																																																																		
<p>THICK. ± 0.05</p> <p>SECTION A-A</p>				<table border="1"> <thead> <tr> <th>PRODUCT #</th> <th>NAME</th> <th>ID</th> <th>OD</th> <th>THICK.</th> <th>MASS (g)</th> </tr> </thead> <tbody> <tr><td>OST05D05</td><td>OPTICAL SPACER RING, 0.5MM THICK, D05</td><td>10.400</td><td>12.600</td><td>0.500</td><td>0.054</td></tr> <tr><td>OST1D05</td><td>OPTICAL SPACER RING, 1MM THICK, D05</td><td>10.400</td><td>12.600</td><td>1.000</td><td>0.105</td></tr> <tr><td>OST2D05</td><td>OPTICAL SPACER RING, 2MM THICK, D05</td><td>10.400</td><td>12.600</td><td>2.000</td><td>0.213</td></tr> <tr><td>OST3D05</td><td>OPTICAL SPACER RING, 3MM THICK, D05</td><td>10.400</td><td>12.600</td><td>3.000</td><td>0.326</td></tr> <tr><td>OST5D05</td><td>OPTICAL SPACER RING, 5MM THICK, D05</td><td>10.400</td><td>12.600</td><td>5.000</td><td>0.551</td></tr> <tr><td>OST05D1</td><td>OPTICAL SPACER RING, 0.5MM THICK, D1</td><td>21.900</td><td>25.300</td><td>0.500</td><td>0.172</td></tr> <tr><td>OST1D1</td><td>OPTICAL SPACER RING, 1MM THICK, D1</td><td>21.900</td><td>25.300</td><td>1.000</td><td>0.335</td></tr> <tr><td>OST2D1</td><td>OPTICAL SPACER RING, 2MM THICK, D1</td><td>21.900</td><td>25.300</td><td>2.000</td><td>0.677</td></tr> <tr><td>OST3D1</td><td>OPTICAL SPACER RING, 3MM THICK, D1</td><td>21.900</td><td>25.300</td><td>3.000</td><td>1.028</td></tr> <tr><td>OST5D1</td><td>OPTICAL SPACER RING, 5MM THICK, D1</td><td>21.900</td><td>25.300</td><td>5.000</td><td>1.706</td></tr> <tr><td>OST05D15</td><td>OPTICAL SPACER RING, 0.5MM THICK, D15</td><td>34.900</td><td>38.000</td><td>0.500</td><td>0.241</td></tr> <tr><td>OST1D15</td><td>OPTICAL SPACER RING, 1MM THICK, D15</td><td>34.900</td><td>38.000</td><td>1.000</td><td>0.471</td></tr> <tr><td>OST2D15</td><td>OPTICAL SPACER RING, 2MM THICK, D15</td><td>34.900</td><td>38.000</td><td>2.000</td><td>0.952</td></tr> <tr><td>OST3D15</td><td>OPTICAL SPACER RING, 3MM THICK, D15</td><td>34.900</td><td>38.000</td><td>3.000</td><td>1.447</td></tr> <tr><td>OST5D15</td><td>OPTICAL SPACER RING, 5MM THICK, D15</td><td>34.900</td><td>38.000</td><td>5.000</td><td>2.401</td></tr> <tr><td>OST05D2</td><td>OPTICAL SPACER RING, 0.5MM THICK, D2</td><td>46.200</td><td>50.700</td><td>0.500</td><td>0.465</td></tr> <tr><td>OST1D2</td><td>OPTICAL SPACER RING, 1MM THICK, D2</td><td>46.200</td><td>50.700</td><td>1.000</td><td>0.913</td></tr> <tr><td>OST2D2</td><td>OPTICAL SPACER RING, 2MM THICK, D2</td><td>46.200</td><td>50.700</td><td>2.000</td><td>1.841</td></tr> <tr><td>OST3D2</td><td>OPTICAL SPACER RING, 3MM THICK, D2</td><td>46.200</td><td>50.700</td><td>3.000</td><td>2.787</td></tr> <tr><td>OST5D2</td><td>OPTICAL SPACER RING, 5MM THICK, D2</td><td>46.200</td><td>50.700</td><td>5.000</td><td>4.630</td></tr> </tbody> </table>				PRODUCT #	NAME	ID	OD	THICK.	MASS (g)	OST05D05	OPTICAL SPACER RING, 0.5MM THICK, D05	10.400	12.600	0.500	0.054	OST1D05	OPTICAL SPACER RING, 1MM THICK, D05	10.400	12.600	1.000	0.105	OST2D05	OPTICAL SPACER RING, 2MM THICK, D05	10.400	12.600	2.000	0.213	OST3D05	OPTICAL SPACER RING, 3MM THICK, D05	10.400	12.600	3.000	0.326	OST5D05	OPTICAL SPACER RING, 5MM THICK, D05	10.400	12.600	5.000	0.551	OST05D1	OPTICAL SPACER RING, 0.5MM THICK, D1	21.900	25.300	0.500	0.172	OST1D1	OPTICAL SPACER RING, 1MM THICK, D1	21.900	25.300	1.000	0.335	OST2D1	OPTICAL SPACER RING, 2MM THICK, D1	21.900	25.300	2.000	0.677	OST3D1	OPTICAL SPACER RING, 3MM THICK, D1	21.900	25.300	3.000	1.028	OST5D1	OPTICAL SPACER RING, 5MM THICK, D1	21.900	25.300	5.000	1.706	OST05D15	OPTICAL SPACER RING, 0.5MM THICK, D15	34.900	38.000	0.500	0.241	OST1D15	OPTICAL SPACER RING, 1MM THICK, D15	34.900	38.000	1.000	0.471	OST2D15	OPTICAL SPACER RING, 2MM THICK, D15	34.900	38.000	2.000	0.952	OST3D15	OPTICAL SPACER RING, 3MM THICK, D15	34.900	38.000	3.000	1.447	OST5D15	OPTICAL SPACER RING, 5MM THICK, D15	34.900	38.000	5.000	2.401	OST05D2	OPTICAL SPACER RING, 0.5MM THICK, D2	46.200	50.700	0.500	0.465	OST1D2	OPTICAL SPACER RING, 1MM THICK, D2	46.200	50.700	1.000	0.913	OST2D2	OPTICAL SPACER RING, 2MM THICK, D2	46.200	50.700	2.000	1.841	OST3D2	OPTICAL SPACER RING, 3MM THICK, D2	46.200	50.700	3.000	2.787	OST5D2	OPTICAL SPACER RING, 5MM THICK, D2	46.200	50.700	5.000	4.630
PRODUCT #	NAME	ID	OD	THICK.	MASS (g)																																																																																																																																
OST05D05	OPTICAL SPACER RING, 0.5MM THICK, D05	10.400	12.600	0.500	0.054																																																																																																																																
OST1D05	OPTICAL SPACER RING, 1MM THICK, D05	10.400	12.600	1.000	0.105																																																																																																																																
OST2D05	OPTICAL SPACER RING, 2MM THICK, D05	10.400	12.600	2.000	0.213																																																																																																																																
OST3D05	OPTICAL SPACER RING, 3MM THICK, D05	10.400	12.600	3.000	0.326																																																																																																																																
OST5D05	OPTICAL SPACER RING, 5MM THICK, D05	10.400	12.600	5.000	0.551																																																																																																																																
OST05D1	OPTICAL SPACER RING, 0.5MM THICK, D1	21.900	25.300	0.500	0.172																																																																																																																																
OST1D1	OPTICAL SPACER RING, 1MM THICK, D1	21.900	25.300	1.000	0.335																																																																																																																																
OST2D1	OPTICAL SPACER RING, 2MM THICK, D1	21.900	25.300	2.000	0.677																																																																																																																																
OST3D1	OPTICAL SPACER RING, 3MM THICK, D1	21.900	25.300	3.000	1.028																																																																																																																																
OST5D1	OPTICAL SPACER RING, 5MM THICK, D1	21.900	25.300	5.000	1.706																																																																																																																																
OST05D15	OPTICAL SPACER RING, 0.5MM THICK, D15	34.900	38.000	0.500	0.241																																																																																																																																
OST1D15	OPTICAL SPACER RING, 1MM THICK, D15	34.900	38.000	1.000	0.471																																																																																																																																
OST2D15	OPTICAL SPACER RING, 2MM THICK, D15	34.900	38.000	2.000	0.952																																																																																																																																
OST3D15	OPTICAL SPACER RING, 3MM THICK, D15	34.900	38.000	3.000	1.447																																																																																																																																
OST5D15	OPTICAL SPACER RING, 5MM THICK, D15	34.900	38.000	5.000	2.401																																																																																																																																
OST05D2	OPTICAL SPACER RING, 0.5MM THICK, D2	46.200	50.700	0.500	0.465																																																																																																																																
OST1D2	OPTICAL SPACER RING, 1MM THICK, D2	46.200	50.700	1.000	0.913																																																																																																																																
OST2D2	OPTICAL SPACER RING, 2MM THICK, D2	46.200	50.700	2.000	1.841																																																																																																																																
OST3D2	OPTICAL SPACER RING, 3MM THICK, D2	46.200	50.700	3.000	2.787																																																																																																																																
OST5D2	OPTICAL SPACER RING, 5MM THICK, D2	46.200	50.700	5.000	4.630																																																																																																																																
<p>UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN MILLIMETERS REMOVE BURRS AND SHARP EDGES TOLERANCES AFTER SURFACE TREATMENT SHOWN: ± 0.1, $\angle \pm 1/4^\circ$ DRAWING IN ACCORDANCE WITH ASME Y14.5 2009</p> <p>MATERIAL: ALUMINUM</p> <p>SURFACE TREATMENT: BLACK ANODIZING</p>		<p>CONTRACT NUMBER/INITIAL PROJECT: 189810</p> <p>APPROVALS: B. LEDUC</p> <p>DATE: 2020-03-30</p> <p>DESIGNED: B. LEDUC</p> <p>DRAWN: B. LEDUC</p> <p>EXAM/CHECK: B. LEDUC</p> <p>2020-04-01</p> <p>2020-04-09</p> <p>2020-04-09</p>		<p>INO NATIONAL OPTICS INSTITUTE 2740, EINSTEIN STREET QUEBEC CITY, QUEBEC CANADA, G1P 4S4</p> <p>T (418) 657-7006 F (418) 657-7009 www.ino.ca</p>		<p>TITLE: OPTICAL SPACER RING</p> <p>FORMAT: A (CAGE CODE) 3AT79</p> <p>DRAWING NUMBER: OST_D_</p> <p>REVISION: A</p> <p>SCALE: 2:1</p> <p>MASS (kg):</p> <p>PAGE: 1/1</p>																																																																																																																															

Optical spacer, special order

Description	Custom-made spacer to hold thicker mirrors, filters, or dichroics in QuickPOZ mirror mount series
Adjustments	n/a
Required tool	n/a
Product notes	Custom spacer available for mount sizes D05, D1, D15, and D2. Custom thickness "THK" shown in drawing section A-A will be determined upon request.

Part no.	Unit Price (\$)
OST_-S	TBA

Drawing OST_-S




Mounting Plates

Breadboard

Description	Mounting plate, 6.35 or 12.7mm thick with M4 x 0.7 threaded holes
Adjustments	n/a
Required tool	n/a
Product notes	Use only 3 of the 8 mounting holes for mounting otherwise the breadboard may warp.

Part no.	Unit Price (\$)
BB064W100L100	TBA
BB064W100L200	TBA
BB064W100L300	TBA
BB127W100L300	TBA
BB127W200L200	TBA
BB127W200L300	TBA
BB127W200L400	TBA
BB127W300L300	TBA
BB127W300L400	TBA

Drawing BB_W_L_




Breadboard disc, gimbal

Description	Mounting disc with rotation adjustment, 6.35mm thick
Adjustments	360°
Required tool	n/a
Product notes	

Part no.	Unit Price (\$)
BBDG1	TBA
BBDG2	TBA

Drawing BBDG_




Gimbal adaptor for mirror mount and fixed mount series

Description	Mounting plate, 6.35mm thick with circular end
Adjustments	360° horizontal travel
Required tool	n/a
Product notes	

Part no.	Unit Price (\$)
BBG150	TBA
BBG200	TBA
BBG250	TBA
BBG300	TBA

Drawing BBG_




Mirror Mounts

Fixed mirror mount, horizontal

Description	Low distortion fixed mirror mount for horizontal beam folding
Adjustments	n/a
Required tool	2 mm Allen wrench
Product notes	<p>Nominal mirror preload is reached at 0.4mm deflection which corresponds to a mirror thickness of:</p> <ul style="list-style-type: none"> • 6 ± 0.1 mm for $\varnothing 12$-12.7mm and $\varnothing 25$-25.4mm mirrors, • 9.525 ± 0.1 mm for $\varnothing 38.1$mm mirrors, • 12 ± 0.1 mm for $\varnothing 50$-50.8mm mirrors. <p>The mount can accommodate other mirror thicknesses by adding optical spacer rings (OST_D_ shims) between the mirror and the blades.</p> <p>If the mirror is just slightly too thick, use ID $\varnothing 1/8''$ x OD $\varnothing 3/16''$ precision shims of the desired thickness from McMaster.</p> <p>For significantly thicker mirrors, use our special OST_-S optical spacer with custom thickness – contact INO for details.</p> <p>*** Warning *** Blades have been designed to be used in stack of 3 with a maximal deflection of 0.5mm at the tip.</p>

Part no.	Unit Price (\$)
MMH254D05	TBA
MMH318D05	TBA
MMH381D05	TBA
MMH254D1	TBA
MMH318D1	TBA
MMH381D1	TBA
MMH254D15	TBA
MMH318D15	TBA
MMH381D15	TBA
MMH318D2	TBA
MMH381D2	TBA

Drawing MMH_D_ (...see next page)



Fixed mirror mount, vertical-bottom

Description	Low distortion fixed mirror mount for vertical beam folding, downwards
Adjustments	n/a
Required tool	2 mm Allen wrench
Product notes	<p>Nominal mirror preload is reached at 0.4mm deflection which corresponds to a mirror thickness of:</p> <ul style="list-style-type: none"> • 6 ± 0.1 mm for $\varnothing 12$-12.7mm and $\varnothing 25$-25.4mm mirrors, • 9.525 ± 0.1 mm for $\varnothing 38.1$mm mirrors, • 12 ± 0.1 mm for $\varnothing 50$-50.8mm mirrors. <p>The mount can accommodate other mirror thicknesses by adding optical spacer rings (OST_D_ shims) between the mirror and the blades.</p> <p>If the mirror is just slightly too thick, use ID $\varnothing 1/8''$ x OD $\varnothing 3/16''$ precision shims of the desired thickness from McMaster.</p> <p>For significantly thicker mirrors, use our special optical spacer OST_-S with custom thickness – contact INO for details.</p> <p>***Warning*** Blades have been designed to be used in stack of 3 with a maximal deflection of 0.5mm at the tip.</p> <p>Optical axis positioning accuracy after the folding may be up to ± 0.1mm RSS.</p>

Part no.	Unit Price (\$)
MMVB254D05	TBA
MMVB254D1	TBA
MMVB254D15	TBA
MMVB318D05	TBA
MMVB318D1	TBA
MMVB318D15	TBA
MMVB381D05	TBA
MMVB381D1	TBA
MMVB381D15	TBA

Drawing MMVB_D_ (...see next page)

4	3	2	1																																																																																																				
PROPRIETARY NOTICE/CONFIDENTIALITY: THIS DOCUMENT AND THE INFORMATION DISCLOSED HEREIN ARE CONFIDENTIAL AND THE PROPERTY OF INSTITUT NATIONAL D'OPTIQUE (NATIONAL OPTICS INSTITUTE) (HEREAFTER "INO"). THEY SHALL NOT BE REPRODUCED OR DISCLOSED IN WHOLE OR IN PART WITHOUT THE EXPRESS WRITTEN AUTHORIZATION OF INO. INO RESERVES ALL PATENT, COPYRIGHT AND OTHER PROPRIETARY RIGHTS TO THIS DOCUMENT, INCLUDING ALL DESIGN, MANUFACTURING, REPRODUCTION, USE AND SALES RIGHTS THEREON. REPRODUCTION, USE AND SALES RIGHTS THEREON, EXCEPT TO THE EXTENT SAID RIGHTS ARE EXPRESSLY GRANTED TO OTHERS.		REVISIONS																																																																																																					
NOTES: 1 - FOR REFERENCE ONLY, NOT FOR MANUFACTURING PURPOSES		REV. DESCRIPTION DATE (YYYY-MM-DD) DRAWN																																																																																																					
<div style="display: flex; justify-content: space-around;"> <div style="text-align: center;"> </div> <div style="text-align: center;"> </div> </div> <div style="text-align: center; margin-top: 10px;"> <p>SECTION A-A</p> </div>																																																																																																							
<table border="1" style="width:100%; border-collapse: collapse; font-size: 8px;"> <thead> <tr> <th>PRODUCT #</th> <th>NAME</th> <th>AXIS HEIGHT</th> <th>INNER DIA.</th> <th>Ø C/A</th> <th>SPACING</th> <th>LENGTH</th> <th>WIDTH</th> <th>HEIGHT</th> <th>MASS (KG)</th> </tr> </thead> <tbody> <tr> <td>MMVB254D05</td> <td>FIXED MIRROR MOUNT, VERTICAL BOTTOM, D05, 25.4MM</td> <td>25.40</td> <td>12.90</td> <td>10.41</td> <td>37.50</td> <td>31.15</td> <td>56.15</td> <td>44.4</td> <td>0.050</td> </tr> <tr> <td>MMVB318D05</td> <td>FIXED MIRROR MOUNT, VERTICAL BOTTOM, D05, 31.8MM</td> <td>31.75</td> <td>12.90</td> <td>10.41</td> <td>37.50</td> <td>31.15</td> <td>56.15</td> <td>50.7</td> <td>0.053</td> </tr> <tr> <td>MMVB381D05</td> <td>FIXED MIRROR MOUNT, VERTICAL BOTTOM, D05, 38.1MM</td> <td>38.10</td> <td>12.90</td> <td>10.41</td> <td>37.50</td> <td>31.15</td> <td>56.15</td> <td>57.1</td> <td>0.056</td> </tr> <tr> <td>MMVB254D1</td> <td>FIXED MIRROR MOUNT, VERTICAL BOTTOM, D1, 25.4MM</td> <td>25.40</td> <td>25.60</td> <td>21.84</td> <td>37.50</td> <td>36.00</td> <td>56.15</td> <td>44.4</td> <td>0.048</td> </tr> <tr> <td>MMVB318D1</td> <td>FIXED MIRROR MOUNT, VERTICAL BOTTOM, D1, 31.8MM</td> <td>31.75</td> <td>25.60</td> <td>21.84</td> <td>37.50</td> <td>36.00</td> <td>56.15</td> <td>50.7</td> <td>0.052</td> </tr> <tr> <td>MMVB381D1</td> <td>FIXED MIRROR MOUNT, VERTICAL BOTTOM, D1, 38.1MM</td> <td>38.10</td> <td>25.60</td> <td>21.84</td> <td>37.50</td> <td>36.00</td> <td>56.15</td> <td>57.1</td> <td>0.057</td> </tr> <tr> <td>MMVB254D15</td> <td>FIXED MIRROR MOUNT, VERTICAL BOTTOM, D15, 25.4MM</td> <td>25.40</td> <td>38.30</td> <td>35.05</td> <td>62.50</td> <td>40.00</td> <td>81.15</td> <td>49.4</td> <td>0.074</td> </tr> <tr> <td>MMVB318D15</td> <td>FIXED MIRROR MOUNT, VERTICAL BOTTOM, D15, 31.8MM</td> <td>31.75</td> <td>38.30</td> <td>35.05</td> <td>62.50</td> <td>40.00</td> <td>81.15</td> <td>55.7</td> <td>0.082</td> </tr> <tr> <td>MMVB381D15</td> <td>FIXED MIRROR MOUNT, VERTICAL BOTTOM, D15, 38.1MM</td> <td>38.10</td> <td>38.30</td> <td>35.05</td> <td>62.50</td> <td>40.00</td> <td>81.15</td> <td>62.1</td> <td>0.089</td> </tr> </tbody> </table>		PRODUCT #	NAME	AXIS HEIGHT	INNER DIA.	Ø C/A	SPACING	LENGTH	WIDTH	HEIGHT	MASS (KG)	MMVB254D05	FIXED MIRROR MOUNT, VERTICAL BOTTOM, D05, 25.4MM	25.40	12.90	10.41	37.50	31.15	56.15	44.4	0.050	MMVB318D05	FIXED MIRROR MOUNT, VERTICAL BOTTOM, D05, 31.8MM	31.75	12.90	10.41	37.50	31.15	56.15	50.7	0.053	MMVB381D05	FIXED MIRROR MOUNT, VERTICAL BOTTOM, D05, 38.1MM	38.10	12.90	10.41	37.50	31.15	56.15	57.1	0.056	MMVB254D1	FIXED MIRROR MOUNT, VERTICAL BOTTOM, D1, 25.4MM	25.40	25.60	21.84	37.50	36.00	56.15	44.4	0.048	MMVB318D1	FIXED MIRROR MOUNT, VERTICAL BOTTOM, D1, 31.8MM	31.75	25.60	21.84	37.50	36.00	56.15	50.7	0.052	MMVB381D1	FIXED MIRROR MOUNT, VERTICAL BOTTOM, D1, 38.1MM	38.10	25.60	21.84	37.50	36.00	56.15	57.1	0.057	MMVB254D15	FIXED MIRROR MOUNT, VERTICAL BOTTOM, D15, 25.4MM	25.40	38.30	35.05	62.50	40.00	81.15	49.4	0.074	MMVB318D15	FIXED MIRROR MOUNT, VERTICAL BOTTOM, D15, 31.8MM	31.75	38.30	35.05	62.50	40.00	81.15	55.7	0.082	MMVB381D15	FIXED MIRROR MOUNT, VERTICAL BOTTOM, D15, 38.1MM	38.10	38.30	35.05	62.50	40.00	81.15	62.1	0.089		
PRODUCT #	NAME	AXIS HEIGHT	INNER DIA.	Ø C/A	SPACING	LENGTH	WIDTH	HEIGHT	MASS (KG)																																																																																														
MMVB254D05	FIXED MIRROR MOUNT, VERTICAL BOTTOM, D05, 25.4MM	25.40	12.90	10.41	37.50	31.15	56.15	44.4	0.050																																																																																														
MMVB318D05	FIXED MIRROR MOUNT, VERTICAL BOTTOM, D05, 31.8MM	31.75	12.90	10.41	37.50	31.15	56.15	50.7	0.053																																																																																														
MMVB381D05	FIXED MIRROR MOUNT, VERTICAL BOTTOM, D05, 38.1MM	38.10	12.90	10.41	37.50	31.15	56.15	57.1	0.056																																																																																														
MMVB254D1	FIXED MIRROR MOUNT, VERTICAL BOTTOM, D1, 25.4MM	25.40	25.60	21.84	37.50	36.00	56.15	44.4	0.048																																																																																														
MMVB318D1	FIXED MIRROR MOUNT, VERTICAL BOTTOM, D1, 31.8MM	31.75	25.60	21.84	37.50	36.00	56.15	50.7	0.052																																																																																														
MMVB381D1	FIXED MIRROR MOUNT, VERTICAL BOTTOM, D1, 38.1MM	38.10	25.60	21.84	37.50	36.00	56.15	57.1	0.057																																																																																														
MMVB254D15	FIXED MIRROR MOUNT, VERTICAL BOTTOM, D15, 25.4MM	25.40	38.30	35.05	62.50	40.00	81.15	49.4	0.074																																																																																														
MMVB318D15	FIXED MIRROR MOUNT, VERTICAL BOTTOM, D15, 31.8MM	31.75	38.30	35.05	62.50	40.00	81.15	55.7	0.082																																																																																														
MMVB381D15	FIXED MIRROR MOUNT, VERTICAL BOTTOM, D15, 38.1MM	38.10	38.30	35.05	62.50	40.00	81.15	62.1	0.089																																																																																														
<div style="display: flex;"> <div style="width: 45%; font-size: 8px;"> <p>UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN MILLIMETERS REMOVE BURRS AND SHARP EDGES TOLERANCES AFTER SURFACE TREATMENT 0.8 / SHOWN ±0.1, <±1/4" DRAWING IN ACCORDANCE WITH ASME Y14.5 2009</p> <p>MATERIAL ALUMINUM</p> <p>SURFACE TREATMENT BLACK ANODIZING</p> </div> <div style="width: 55%; font-size: 8px;"> <p>CONTRACT NUMBER/INITIAL PROJECT 189810</p> <p>APPROVALS DATE (YYYY-MM-DD)</p> <p>DESIGNED M.GRENIER 2020-02-05</p> <p>DRAWN M.GRENIER 2020-02-24</p> <p>EXAM/CHECK <i>SG</i> 2020-04-30</p> <p>VÉRIF./VÉRIF. <i>MD</i> 2020-04-30</p> </div> </div>		<div style="display: flex; justify-content: space-between;"> <div style="width: 40%; font-size: 8px;"> <p>INO ALL RIGHTS RESERVED ©</p> <p>NATIONAL OPTICS INSTITUTE 2740, EINSTEIN STREET QUEBEC CITY, QUEBEC CANADA, G1P 4S4</p> </div> <div style="width: 55%; font-size: 8px;"> <p>T (418) 657-7006 F (418) 657-7009 www.ino.ca</p> </div> </div> <p>TITLE FIXED MIRROR MOUNT, VERTICAL-BOTTOM</p> <div style="display: flex; justify-content: space-between; font-size: 8px;"> <div> <p>FORMAT DAI (CAGE CODE) A 3AT79</p> <p>SCALE 1:2</p> </div> <div> <p>DRAWING NUMBER MMVB_D_</p> <p>MASS (kg) - PAGE 1/1</p> </div> <div> <p>REV. A</p> </div> </div>																																																																																																					
4	3	2	1																																																																																																				

Adjustable mirror mount, horizontal

Description	Low distortion adjustable mirror mount with $\pm 2^\circ$ tip-tilt for horizontal beam folding. Left-hand and right-hand versions.
Adjustments	D05 & D1: travel tip/tilt $\pm 2^\circ$; Z ± 1 mm, 8.4 mrad/rev. D15: travel tip/tilt $\pm 2^\circ$; Z ± 1 mm, 6.6 mrad/rev. D2: travel tip/tilt $\pm 2^\circ$; Z ± 1 mm, 5.3 mrad/rev. Resolution after locking: 0.010-0.015 mrad.
Required tool	2mm Allen wrench
Product notes	<p>A 2 or 2.5 mm gage pin can be used for nominal gap axial control position.</p> <p>Belleville spring stack nominal position is reached at one turn of M4 screw loosening with respect to fully compress position.</p> <p>For locking, fully compress Belleville spring stack and then loosen for $\frac{1}{4}$ turn of M4 screw.</p> <p>Nominal mirror preload is reached at 0.4mm deflection which corresponds to a mirror thickness of:</p> <ul style="list-style-type: none"> • 6 ± 0.1 mm for $\varnothing 12$-12.7mm and $\varnothing 25$-25.4mm mirrors, • 9.525 ± 0.1 mm for $\varnothing 38.1$mm mirrors, • 12 ± 0.1 mm for $\varnothing 50$-50.8mm mirrors. <p>The mount can accommodate other mirror thicknesses by adding optical spacer rings (OST_D_ shims) between the mirror and the blades.</p> <p>If the mirror is just slightly too thick, use ID $\varnothing 1/8''$ x OD $\varnothing 3/16''$ precision shims of the desired thickness from McMaster.</p> <p>For significantly thicker mirrors, use our special optical spacer OST_-S with custom thickness – contact INO for details.</p> <p>***Warning*** Blades have been designed to be used in stack of 3 with a maximal deflection of 0.5mm at the tip.</p>

Part no.	Unit Price (\$)
MMALH318D05	TBA
MMALH318D1	TBA
MMALH381D05	TBA
MMALH381D1	TBA
MMALH381D15	TBA
MMALH381D2	TBA
MMARH318D05	TBA
MMARH318D1	TBA
MMARH381D05	TBA
MMARH381D1	TBA
MMARH381D15	TBA
MMARH381D2	TBA

Drawing MMA_H_D_ (...see next page)

Released A 20. MMA H D


Generic Mounts

XY adjustable mount

Description	X-Y adjustable mount with axial support and QC threads
Adjustments	Refer to TLXY1
Required tool	2mm Allen wrench
Product notes	Compatible with QuickPOZ TLXY1 adjustment tool, QuickPOZ Tube QC series, commercial threaded and unthreaded accessories.

Part no.	Unit Price (\$)
MA254QC05	TBA
MA318QC05	TBA
MA381QC05	TBA
MA254QC1	TBA
MA318QC1	TBA
MA381QC1	TBA

Drawing MA_QC_




Through-hole XY adjustable mount

Description	X-Y adjustable mount with axial support and through-hole
Adjustments	Refer to TLXY1
Required tool	2mm Allen wrench
Product notes	Compatible with QuickPOZ TLXY1 adjustment tool, QuickPOZ Tube QC thread series, commercial threaded and unthreaded accessories.

Part no.	Unit Price (\$)
MA254TH05	TBA
MA318TH05	TBA
MA381TH05	TBA
MA254TH1	TBA
MA318TH1	TBA
MA381TH1	TBA

Drawing MA_TH_




XYZ adjustable mount

Description	X-Y-Z adjustable mount with QC threads
Adjustments	Refer to TLXY1 for X & Y adjustment. For axial adjustment: 7mm travel, 2µm axial resolution (1° of adjuster rotation), Knurled locking ring.
Required tool	2mm Allen wrench
Product notes	Compatible with QuickPOZ TLXY1 adjustment tool, QuickPOZ Tube QC series, commercial threaded and unthreaded accessories.

Part no.	Unit Price (\$)
MA254TAQC1	TBA
MA318TAQC1	TBA
MA381TAQC1	TBA

Drawing MA_TAQC_




Fixed mount with flange

Description	Fixed mount with axial support and QC threads
Adjustments	n/a
Required tool	n/a
Product notes	Compatible with QuickPOZ QC threads as well as commercial threads

Part no.	Unit Price (\$)
MF254QC05	TBA
MF318QC05	TBA
MF381QC05	TBA
MF254QC1	TBA
MF318QC1	TBA
MF381QC1	TBA
MF254QC30	TBA
MF318QC30	TBA
MF381QC30	TBA
MF318QC2	TBA
MF381QC2	TBA

Drawing MF_QC_




Fixed mount

Description	Fixed mount with QC threads
Adjustments	n/a
Required tool	n/a
Product notes	Compatible with QuickPOZ QC threads as well as commercial threads

Part no.	Unit Price (\$)
M254QC05	TBA
M318QC05	TBA
M381QC05	TBA
M254QC1	TBA
M318QC1	TBA
M381QC1	TBA
M254QC30	TBA
M318QC30	TBA
M381QC30	TBA
M318QC2	TBA
M381QC2	TBA
M254RMS	TBA
M318RMS	TBA
M381RMS	TBA

Drawing M_QC_



Tubes


Lens tube, male-female

Description	Autocentered tube with QC threads
Adjustments	n/a
Required tool	n/a
Product notes	Compatible with QuickPOZ QC threads as well as commercial threads

Part no.	Unit Price (\$)
TMF064QC05	TBA
TMF127QC05	TBA
TMF254QC05	TBA
TMF381QC05	TBA
TMF508QC05	TBA
TMF064QC1	TBA
TMF127QC1	TBA
TMF254QC1	TBA
TMF381QC1	TBA
TMF508QC1	TBA

Part no.	Unit Price (\$)
TMF064QC30	TBA
TMF127QC30	TBA
TMF254QC30	TBA
TMF381QC30	TBA
TMF508QC30	TBA
TMF064QC2	TBA
TMF127QC2	TBA
TMF254QC2	TBA
TMF381QC2	TBA
TMF508QC2	TBA

Drawing TMF_QC_




Lens tube, male-female, adapted for MA_TH_

Description	Autocentered tube with QC threads, adjusted for use with MA_TH_ mount
Adjustments	n/a
Required tool	n/a
Product notes	Compatible with QuickPOZ QC threads as well as commercial threads

Part no.	Unit Price (\$)
TMFR244QC05	TBA
TMFR244QC1	TBA

Drawing TMFR_QC_




Lens tube with axial adjustment

Description	Tube with axial adjustment
Adjustments	7 mm travel, 2µm axial resolution (1° of adjuster rotation), Knurled locking ring.
Required tool	n/a
Product notes	Compatible with QuickPOZ QC threads as well as commercial threads

Part no.	Unit Price (\$)
TAQC05	TBA
TAQC1	TBA
TAQC30	TBA
TAQC2	TBA

Drawing TAQC_




Rotation Mounts

Rotation mount

Description	Rotation mount
Adjustments	Travel: 360° endless; resolution $\pm 1^\circ$, self locking.
Required tool	Pin $\varnothing 1$ mm
Product notes	Compatible with QuickPOZ QC threads as well as commercial threads

Part no.	Unit Price (\$)
RMQC05M	TBA
RMQC05MF	TBA
RMQC1M	TBA
RMQC1MF	TBA

Drawing RMQC_




Translation Mounts

Translation mount, 1.5mm, 62.5MM

Description	Translation flexure, $\pm 1.5\text{mm}$, 62.5MM
Adjustments	$\pm 1.5\text{mm}$; 300 μm /revolution, resolution of $\pm 0.8\mu\text{m}$
Required tool	2mm Allen wrench
Product notes	Compatible with QuickPOZ MA_QC05 & MA_QC1 mounts. The use of Belleville springs is recommended to spring load QuickPOZ MA_QC_mount during adjustment.

Part no.	Unit Price (\$)
MZ1	TBA

Drawing MZ1





INO
2740 Einstein Street
Quebec City, Quebec
G1P 4S4 Canada

www.ino.ca

418-657-7006 / 1-866-657-7406
info@ino.ca